

Automated Scan Downloading and Project Organisation

23042021


Overview:

- Initialise project in pre-determined BIDS folder structure
- Download scans from WBIC, organise with Scan ID (potential for automatic backups)
- Copy user-specified raw scans to pre-processing folder
- Users only need to specify options in a configuration file – no coding required
- Object-orientated approach means easy to duplicate for new projects and customise
- Full logging and testing

Current Features:

- Download scans from WBIC to HPC then from HPC to Hivemind
- Deletes data from HPC to avoid data buildup
- Reorganise downloaded data to ABL format
- Copy raw data to preprocessing with flexible user-specified search strings.
Sessions are automatically added and text files including all important subject IDs, scan date, time etc. written to session folder
- Optional flags for skipping runs (e.g. if a run fails)
- Every command is logged (e.g copy file, move file) for convenient checking
- All user input is tested to ensure no input errors
- Scan / Session IDs are tested against scan time / date to ensure order is corrected
- Maintains BIDS organisation throughout project

Project / BIDS structure



/mnt/zfs/PROJECT_NAME/ – docs

– logs

data

– pre-processing

– sub-00X

– ses-00X

– mrs

– raw

– sub-00X_task_XXX_run-X_slaser

– .dcms

– func ...


– anat ...

– mpm ...

– ses-00x_info.txt (WBIC ID, ABL ID, sub, ses ID, scan zk ID, scan date, scan time)

– raw_scans - zk.....001

- zk-....002



Name	Date modified	Type	Size
DATA_0039_00001_0000.dcm	23/04/2021 12:23	Adobe Photoshop...	137 KB
DATA_0039_00002_0000.dcm	23/04/2021 12:23	Adobe Photoshop...	137 KB
DATA_0039_00003_0000.dcm	23/04/2021 12:23	Adobe Photoshop...	137 KB
DATA_0039_00004_0000.dcm	23/04/2021 12:23	Adobe Photoshop...	137 KB
DATA_0039_00005_0000.dcm	23/04/2021 12:23	Adobe Photoshop...	137 KB
DATA_0039_00006_0000.dcm	23/04/2021 12:23	Adobe Photoshop...	137 KB

Code Organisation


/git-repos/fmri_scan_fetch/

- backend
 - all backend code
- project_configs.py (only thing user needs to change)


ProjectMaster superclass

Project subclass that inherits ProjectMaster methods and overwrites config attributes.

This means full generalisation across projects and easy customisation by method subclassing



Name	Date modified
__pycache__	20/04/2021 23:32
backend	23/04/2021 12:47
__init__.py	20/04/2021 13:09
project_configs.py	23/04/2021 12:10
project_configs.pyc	20/04/2021 20:21



Name	Date modified
__pycache__	23/04/2021 12:47
__init__.py	20/04/2021 13:09
__init__.pyc	20/04/2021 20:21
auto_preprocess_wbic_scans.py	23/04/2021 11:11
project_master.py	23/04/2021 11:11
project_master.pyc	20/04/2021 20:21

Config File – Project Details

<pre># Project Settings ----- self.base_path = "/mnt/zfs/7T_TUNING" self.scripts_path = "/mnt/Data/jjz33/hpc/backend" self.project_code = "P00461" self.account = "jjz33" self.scanner_format = ".dcm" self.mrs_scan_info = {"slaser": {"search_str": "*_slaser_WS128_OCC_w_Pad_LongTE", "task_name": "ori"} }</pre>	<p>Project details</p> <p># Base directory for the folder</p> <p>#</p> <p># WBIC project code for the project</p> <p># User ID. This should be the same for the Hivemind / HPC / WBIC</p> <p># search_str used with glob.glob, see docs for it's format options (its different to regexp)</p>
<pre>self.func_scan_info = {"bold_ptx": {"search_str": "*_rs1h_ep3d_vaso_ma4a_TReff4s_30sl_Run?_E00_M", "task_name": "ori"}, "bold_old": {"search_str": "*_rs1h_ep3d_vaso_ma4a_TReff4s_30sl_R?_E00_M", "task_name": "ori"} }</pre>	<p>< Flexible string searching for files to copy from raw scans to preprocesing</p> <p>< handles different tasks / runs / scan types</p>
<pre>self.anat_scan_info = {"mp2rage": {"search_str": "*mp2rage_sag_p2_0.65mm_UNI_Images", "task_name": "ori"} } self.mpm_scan_info = None</pre>	<p>Input expected number of volumes for testing</p> <pre>self.num_expected_func_files = 292 self.num_expected_anat_files = 240 self.num_expected_mrs_files = 136 self.num_expected_mpm_files = None</pre>

Config File – Participants

```
# Participants -----

self._participant_log = {

    "31970": {"sub_id": "sub-001", All sub IDs in one place
        "lab_id": "3470",
        "scans": {"scan_1": {"date": "20210420", Supports multiple scans per participant
            "zk_id": "zk21w7_006",
            "time_start": "09:30",
            "flags": ["ignore_mrs_1"]}, Indicate bad runs
        },
    },

    "31974": {"sub_id": "sub-002",
        "lab_id": "3471",
        "scans": {"scan_1": {"date": "20210420",
            "zk_id": "zk21w7_007",
            "time_start": "12:45",
            "flags": ["ignore_mrs_1"],
        },
    },

    "31879": {"sub_id": "sub-003",
        "lab_id": "3472",
        "scans": {"scan_1": {"date": "20210421",
            "zk_id": "zk21w7_008",
            "time_start": "14:00"}
        },
    },
}
```

Logs all SSH connection information for easy troubleshooting

```

Logger initialised. Pulling scans from hpc...
starting thread (client mode): 0x176c6ba8
Local version/idstring: SSH-2.0-paramiko 2.7.2
Remote version/idstring: SSH-2.0-OpenSSH 7.4
Connected (version 2.0, client OpenSSH 7.4)
kex algos:['curve25519-sha256', 'curve25519-sha256@libssh.org', 'ecdh-sha2-nistp256', 'ecdh-sha2-ssh.com', 'umac-128-etm@openssh.com', 'hmac-sha2-256-etm@openssh.com', 'hmac-sha2-512-etm@openssh.com']
Kex agreed: curve25519-sha256@libssh.org
HostKey agreed: ssh-ed25519
Cipher agreed: aes128-ctr
MAC agreed: hmac-sha2-256
Compression agreed: none
kex engine KexCurve25519 specified hash_algo <built-in function openssl_sha256>
Switch to new keys ...
Adding ssh-ed25519 host key for login.hpc.cam.ac.uk: b'ebe3a1f06468cf9c63da8db2ee1583'
Trying SSH key b'6767ddaab4c9f4b757b6f4fa6ebb2c6a'
userauth is OK
Auth banner: b'\n          <><><><><><><><><><><><><><><><><><>\n'
Authentication continues...
Methods: ['keyboard-interactive', 'hostbased']
userauth is OK
Authentication (keyboard-interactive) successful!
[chan 0] Max packet in: 32768 bytes
Received global request "hostkeys-00@openssh.com"
Rejecting "hostkeys-00@openssh.com" global request from server.
[chan 0] Max packet out: 32768 bytes
Sesch channel 0 opened.
[chan 0] Sesch channel 0 request ok
[chan 0] EOF received (0)
[chan 0] EOF sent (0)
Dropping user packet because connection is dead.


Fri Apr 23 00:00:00 2021 pulled scans from wbic to hpc -----
project: P00461, wbic_id 31970, date: jjz33
Patient ID Date Study Description Acc. Number
1 31970 20210420 WBIC Protocols Ready U-ID51976

a: 1 AAHead_Scout_32ch-head-c
b: 2 AAHead_Scout_32ch-head-c
c: 3 AAHead_Scout_32ch-head-c
d: 4 AAHead_Scout_32ch-head-c
e: 5 bl_mapping_2mm
f: 6 bl_mapping_2mm
g: 7 mp2rage_sag_p2_0.65mm_IN
h: 8 mp2rage_sag_p2_0.65mm_IN
i: 9 mp2rage_sag_p2_0.65mm_UN
j: 10 mp2rage_sag_p2_0.65mm_UN
k: 11 t2_tse_tra_fast_for_repo
l: 12 gre_b0map_more_slices
m: 13 gre_b0map_more_slices
n: 14 rs1h_ep3d_vaso_ma4a_TRef
o: 15 AAHead_Scout_32ch-head-c
p: 16 rs1h_ep3d_vaso_ma4a_TRef
q: 17 rs1h_ep3d_vaso_ma4a_TRef

```

Every file operation is logged, and tests reported

```
Fri Apr 23 00:00:00 2021 Processed a new session -----
sub: sub-003,
existing ses files: []
new_ses_id: ses-001

Fri Apr 23 00:00:00 2021 copying raw data to preprocessing folder -----
data type: mrs, zk_id: zk21w7_008, sub_id: sub-003

copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_039_slaser_WS128_OCC_w_Pad_LongTE/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/mrs/raw/sub-003_task-ori_run-1_slaser
scan has the expected number of files: 136

Fri Apr 23 00:00:00 2021 copying raw data to preprocessing folder -----
data type: func, zk_id: zk21w7_008, sub_id: sub-003

copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_015_rslh_ep3d_vaso_ma4a_TReff4s_30s1_Run1_E00_M/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-1_bold_ptx
WARNING: Only 28 scans in /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-1_bold_ptx but expecting 292

copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_017_rslh_ep3d_vaso_ma4a_TReff4s_30s1_Run1_E00_M/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-2_bold_ptx
scan has the expected number of files: 292

copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_018_rslh_ep3d_vaso_ma4a_TReff4s_30s1_Run2_E00_M/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-3_bold_ptx
scan has the expected number of files: 292

copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_019_rslh_ep3d_vaso_ma4a_TReff4s_30s1_Run3_E00_M/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-4_bold_ptx
scan has the expected number of files: 292

copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_020_rslh_ep3d_vaso_ma4a_TReff4s_30s1_Run4_E00_M/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-5_bold_ptx
scan has the expected number of files: 292


copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_021_rslh_ep3d_vaso_ma4a_TReff4s_30s1_Run5_E00_M/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-6_bold_ptx
scan has the expected number of files: 292

copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_043_rslh_ep3d_vaso_ma4a_TReff4s_30s1_Run5_E00_M/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-7_bold_ptx
WARNING: Only 10 scans in /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-7_bold_ptx but expecting 292
```

Warnings when something goes wrong

```
copied from: /mnt/zfs/7T_TUNING/data/raw_scans/zk21w7_008/Series_015_rslh_ep3d_vaso_ma4a_TReff4s_30sl_Run1_E00_M/*
copied to: /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-1_bold_ptx
WARNING: Only 28 scans in /mnt/zfs/7T_TUNING/data/preprocessing/sub-003/ses-001/func/raw/sub-003_task-ori_run-1_bold_ptx but expecting 292
```


Folder Organisation and Session info file



BIDS folder structure is imposed

< all session details collated in single file per session

< scan date / times used to test subject order



BIDS file structure

Next: Automated early pre-processing steps in AFNI, FSL, SPM